

## ID'SPASSOR

THE SOLUTION TO MANAGE THE VOLUME OF THE PROSTHETIC TOOTH.

### CONTENTS :

6 pointing drills.  
Drilling part length :  
7 mm and Ø2 mm.

6 spacing guides.  
From Ø4,5 mm to Ø12 mm.



### ID-SPASSOR PROTOCOL

#### STEP 1

Use the pointing drill adapted to the mesio distal space of the tooth and stick it to the side of the adjacent tooth, then proceed to drilling.



#### STEP 2

Remove the drill and replace it by the spacing guide of the same color.



**The drilling guides** are intended to create a predrilling position for an implant setting, perfectly positioned in the middle of the future prosthetic tooth.

**The spacing guide** materializes the perfect volume of the prosthetic tooth.

1. Drill with the pointing drill

2. Remove the pointing drill

3. Insert the spacing guide

